## Imaging in cancer immunology: Phenotyping of multiple immune cell subsets in situ in FFPE tissue sections

Maciej P Zerkowski

PerkinElmer, USA, E-mail: Maciej.Zerkowski@PERKINELMER.COM

## Abstract

There has been a rapid grown in the field of tumor immunobiology in recent years as a result of recent successes in cancer immunotherapies and it is becoming clear that immune cells play many sometimes conflicting roles in the tumor microenvironment. However, obtaining phenotypic information about the various immune cells that play these roles in and around the tumor has been a challenge. Existing methods can either deliver phenotypic information on homogenous samples (e.g., flow cytometry PCR) or morphologic information on single or immunomarkers (standard IHC). We present here a methodology for delivering quantitative per-cell marker expression and phenotyping, analogous to that obtained from flow cytometry but from cells imaged in situ in FFPE tissue sections. This methodology combines the sequential multi-marker labeling of up to 6 antigens using antibodies all of the same species in a single section; automated multispectral imaging (MSI) to remove the typically problematic FFPE tissue auto fluorescence and correct cross-talk between fluorescent channels and an automated image analysis that can quantitate the per-cell marker expression, determine the cellular phenotype, count these cells separately in the tumor compartment and in the stroma and provide high-resolution images of their distributions.

## **Biography**

It have more than 15 years of experience in oil and gas industry (incl. Subsea, off-shore and on-shore) as well as other industries in issue such as fracture mechanics, forensic engineering, integrity, corrosion management, FEED, design and materials selection, both in Australia and overseas. With a background in Metallurgy and Materials Science and Engineering, I have experience in weld inspection, mechanical failure assessment (HISC, SCC, SSCC, MICSCC) microbial corrosion treatment, biocide/inhibitor selection and application and employment of CRA, CP design, coating evaluation for both sour and sweet corrosion conditions.

We present here several examples of this new methodology in breast, lung and head and neck cancers. Each application example will show 6-plex multiplexed staining, per-cell quantitation of each marker and multi-marker cellular phenotyping from multispectral images of standard clinical biopsy sections as well as methods to explore the spatial distributions of the phenotyped cells in and around the tumor.

This work is partly presented at

5th International Conference on Pathologyon May 9-11, 2016 held in Chicago IL, USA.