Development and clinical validation of antioxidants based dried blood spots assay - Status as early biochemical marker for rare genetic disorders in Neonates.

Shyama Subramaniam¹, Samu Subramaniam^{2*}, Subhasree. C.R.³, Arasar Seeralar. A.T.⁴

¹Professor, Consultant, Lab services, Apollo Hospitals, Chennai, Tamil Nadu, India

²Lab Director and Head of the Department, Department of Biochemistry, Regenix Super Specialty Laboratories Private Ltd, Chennai, India

³Research Scholar, Department of Biochemistry, Regenix Super Specialty Laboratories Private Ltd, affiliated to University of Madras, Chennai, Tamil Nadu, India

⁴Professor of Neonatology, Director and Superintendent, Institute of Child Health and Hospital for Children, Madras Medical College, Egmore, Chennai, Tamilnadu, India

Abstract

Introduction: The present study is to investigate the total antioxidant activity (TAC) in suspected neonates of inborn errors of metabolism with high risk, moderate risk and low risk. Dried blood spots have potential use in remote health applications for individual and population diagnosis and can enable epidemiological surveillance for known and unknown diseases. Dried blood spot sampling is most common method for collection, storing, transporting and analyzing a variety of human body fluids.

Materials and methods: After performing, while newborn screening is early period, the balance residual samples were used for specific secondary research studies with patient consent. Here we are using the secondary studies of DBS for total antioxidant activity as one of the biomarkers to find out the disease condition. TAC was evaluated in 449 risk subjects and 531 control neonates including high, moderate, low risk and healthy individuals. The risk variables were included birth weight, APGAR score, gestational age and previous complications of pregnancy. DBS -TAC method, is to elute sample by using phosphate buffered saline. After sample preparation, TAC by using FRAP assay (ferric reducing antioxidant power assay) were estimated by among four groups with risk indicators.

Results: Statistical comparisons and correlations at 5% level of significance were determined. The mean TAC concentration was significantly elevated in high risk group. The mean TAC for low, moderate and high-risk group (Mean \pm SD values were 954.8 \pm 132.2 µmol/L, 700.4 \pm 100.6 µmol/L and 510.4 \pm 110 µmol/L respectively) patients were significantly reduced compared with control group (Mean \pm S.D. valued was 1017.3 \pm 122.2 µmol/L) (p < 0.05). Likely the SOD, Vit C, Vit E and MDA also were estimated.

Discussion: DBS-TAC can be used as early biochemical marker for oxidative stress in high risk for IEM (Inborn errors of metabolism) which may result in reduced tissue damage by free radicals and help to monitor and optimize antioxidant therapy in such high risk neonates.

Conclusion: FRAP assay is presented as a novel method of assessing total antioxidant capacity and is considered as a useful indicator of the body's antioxidant status to counteract the oxidative damage due to ROS. The advantage of the FRAP assay is in being fast, easy to handle, with highly reproducible results. A significantly lower DBS TAC concentration in the IEM high risk patients compared to controls reflects a lower total antioxidant capacity. Thus, TAC may be useful as an early marker of oxidative stress to monitor and optimize antioxidant therapy as an adjunct in the management of IEM patients.

Keywords: Antioxidant therapy, Lipd peroxidation, Metabolic disorders, Newborn Screening, Superdioxide mutase, Vitamin C, Vitamin E.

^{*}Correspondence to: Samu Subramaniam, Departments of Biochemistry, Lab Director and Head of the Department, Regenix Super Specialty Laboratories Private Ltd, Chennai, India, E-mail: Subramaniam@regenix.in

Received: 17-May-2022, Manuscript No. AACBC-22-64084; Editor assigned: 20-May-2022, PreQC No. AACBC-22-64084(PQ); Reviewed: 16-Jun-2022, QC No. AACBC-22-64084; Published: 30-Jun-2022, DOI:10.35841/aacbc-6.3.111

Introduction

Dried blood sampling and newborn screening aims at the earliest possible recognition of disorders to prevent the most serious consequences by timely intervention screening is not a confirmatory diagnosis and requires further investigations. New-born screening is a simple blood test to identify neonates who may have one of 29 primary treatable rare inherited disorders. A new-born baby can look healthy but one of these rare disorders. The most common disorders are congenital hypothyroidism (CH), cystic fibrosis (CFTR), medium chain acyl coA deficiency (MCAD), maple syrup urine disease (MSUD), congenital adrenal hyperplasia (CAH), biotinidase deficiency (BIOT), phenylketonuria (PKU), galactosemia (GALT) and glucose-6-phosphatase deficiency (G-6-PD).

These diseases will be included in the first case, because the annual birth rate is 21.76 births /1000 population and in Delhi alone, nearly 900 births takes place every day; considering this figure there may be one or two babies born Delhi alone with a metabolic defect each day [1]. One in every 1,000 newborns has inborn errors of metabolism (IEM) or other disorders that can be detected at birth by screening (**Table 1**). Inborn errors of metabolic disorders are rare genetic (inherited) disorders in which the body cannot properly turn food into energy. The disorders are usually caused by defects in specific proteins (enzymes) that help break down (metabolize) parts of food [2].

In the present study, we investigated total antioxidant activity as early biochemical marker of oxidative damage in dried blood spots samples from IEM subjects with high, moderate and low risk to measure of lipid peroxidation and super dioxide mutase, a measure of protein oxidation. We also determined DBS- total antioxidant activity in these samples. Subsets of IEM patients showed significantly higher levels of the damage markers with total antioxidant activity.

Study Design

The study has been approved by the institutional human ethical committee approval from Madras Medical College, Chennai: Approval No. 32042018 and obtained consent form from parents for the participation of this research study. Under aseptic conditions with precautions about dried blood spots sampling was collected from suspected IEM neonates and from healthy controls. All neonates aged (0-28) days after birth, normal newborns and including low birth weight neonates and that admitted to Neonatal Intensive Care Unit (NICU), Institute of obstetrics and Gynaecology and Govt. hospital for Women and Children, will be systematically enrolled into the study. The neonatal heel prick or Guthrie test is a screening test done on newborns. It consists of making a pinprick puncture in one heel of the newborn and soaking the blood into pre-printed collection cards known as Guthrie cards. Heel prick capillary blood collected on S&S 903 filter paper, air dried and transported along with ice pack in cold storage.

Methods and Materials

The study was conducted at Department of Biochemistry, Regenix Super speciality laboratories pvt ltd, Chennai. The study included 980 neonates aged under 28 days of life of either sex (201 male babies with mean age of 11 days, 223 female babies with mean age of 9 days and 25 ambigous genitalia babies with mean age of 6 days) were suspected as inherited metabolic disorders confirmed by specific assays like BIORAD-Quantase neonatal screening, aminoacid and acyl carnitine profiling by Tandem mass spectrometry, before the commencement of therapy and 531 neonates, age and sex matched healthy controls. The demographics of this study population shows risk factors associated with maternal complications type of delivery, familial conditions, mother's history, drugs usage, gestational age, number of births, and number of previous siblings death (Figure 1). IEM suspected neonates were further classified according to the ACMG recommended panel criteria. All subjects underwent an initial screening that included a detailed history, signs, symptoms, clinical and laboratory findings of neonates were eligible if they had symptomatic conditions. Exclusion criteria were as follows as: signs of severe infection, severe vomiting, and congenital anomalies.

Modified Method of Total Antioxidant activity -Ferric reducing power assay by using dried blood spots

Elution Tactics

Phosphate Buffered with Saline (10X): 800 ml of distilled water; 8 g of Nacl; 200 mg of Kcl; 1.44 g of Na₂HPO₄; 240 mg of KH₂PO₄, pH-7.4 and add distilled water upto 1L. For sample preparation, simply punch out single 3 mm disc with dried blood spots in a 1.5 ml eppendrof, Add 200 ul of phosphate buffered saline then centrifuge with 12,000 rpm spin column for 2 to 5 mintues. The eluted sample is ready to use. The standard method for the total antioxidant activity was estimated by FRAP method by Benzie.F.F. and J.J.Strain. Total antioxidant capacity was measured by FRAP assay according to the method of Benzie.F.F. and J.J.Strain [3]. At low pH, when a ferric tripyridyltriazine (Fe III-TPTZ)

Table 1. Burden of Genetic and Metabolic Diseases in India: Recent Studies (2015).

Disorder	Incidence	Births/year
Congenital Malformation	1:50	678,000
Down's Syndrome	1:800	34,000
Metabolic Disorders	1:1200	22,500
Beta-Thalassemia and Sickle cell Disease	1:1700	16,700
Congenital Hypothyroidism	1:2500	10,900
Duchenne muscular dystrophy	1:10000	2,700
Spinal Muscular atrophy	1:10000	2,700

complex is reduced to the ferrous (Fe II) form, an intense blue colour with an absorption maximum at 593 nm develops. Antioxidant power convert ferric to ferrous tripyritriazine complex. Ferric reducing ability of plasma (FRAP) values are obtained by comparing the absorbance change at 593 nm in mixture(test), with those containing ferrous ion in known concentration (**Table 2**).

Procedure

Acetate buffer (pH 3.6): 3.1 gm sodium acetate was dissolved in 500 ml distilled water and to this 16 ml of glacial acetic acid was added. pH of the solution was adjusted up to 3.6 by 1% Na₂Co₃ and finally diluted with distilled water to 1000 ml. Ferric Chloride (20M/I): 324 mg. of Fecl₃.6H₂O was dissolved in distilled water and diluted to 100 ml. Ferrozine (10 mM): 49 mg Ferrozine was dissolved in 10 ml of 0.04 M HCL. Standard (1mmol/L): 278 mg Iron(II) sulfate heptahydrate (FeSO₄.7H₂O) was dissolved in 1000 ml distilled water.

Like FRAP assay, estimation of lipid peroxidation by MDA, Vit E, Vit C and SOD enzyme levels were carried out on DBS samples were measured by UV-Spectrophotometry (Chemito). Vitamin E was estimated by Baker and Frank method [4] and vitamin C by 2, 4 – DNPH method [5]. SOD enzyme activity was determined with the direct spectrophotometric method employing KO2 as previously described by Marklund, with some modifications [6,7]. Lipid peroxidation was measured by MDA estimation according to the colorimetric method of Satoh.k [8]. Lipoproteins are precipitated from the specimen by adding trichloroacetic acid. 0.05 M sulphuric acid and 0.67% thiobarbituric acid (TBA) in 2 M sodium sulphate are added to this precipitate and the coupling of lipid peroxide with TBA is carried out by heating in a boiling water bath for 30 minutes. The resulting chromogen is extracted in n-butanol, which is measured at 530 nm.

Statistical Analysis

Each result was expressed as mean \pm standard deviation. The statistical significance of the data was determined by Student's t-test and one way ANOVA test at 5% level of significance (**Table 3**). Pearson's correlation coefficient was determined at 5% level of significance. Statistical analysis was done using SPSS software version 21.0.

Results

There was a significantly increased levels of MDA in the high risk patients comparing to moderate, low risk and control groups (p<0.05). They were significantly decrease levels of Vitamin E, Vitamin C, SOD enzyme activity levels, along with the TAC in the positive patients (significantly decreased in IEM high risk positive patients compared with the control groups(p<0.05). Statistically significant negative correlations found between TAC & MDA in high and moderate risk




Figure 1. Demographics features of the study population: Maternal complications with Low/Moderate/ High risk neonates.

Table 2.	New Assay method for	·DBS-TAC (dried	blood spots- total an	tioxidant activity by FRAP	assay).
----------	----------------------	-----------------	-----------------------	----------------------------	---------

500 ul 50 ul 50 ul 100ul	550 ul 50 ul	500 ul 50 uL 50 ul	500 ul 50 ul 50 ul
50 ul	50 ul		
		50 ul	50 ul
100ul			
10041			
		100 ul	
	100 ul		
Mixe	ed well and incubated at 37 C for 3	30 min	
	100 ul		
	 Mixe	100 ul Mixed well and incubated at 37 C for 3 100 ul	100 ul Mixed well and incubated at 37 C for 30 min

Antioxidant status	High Risk	Moderate Risk	Low Risk	Controls	P value
TAC(FRAP) (µ mol/L)	510.4 ± 110.1	700.4 ± 100.6	954.8 ± 132.2	1017.3 ± 122.2	p<0.05
SOD (units/ml)	3.4 ± 0.4	4.1 ± 0.4	5.7 ± 1.1	7.7 ± 0.9	p<0.05
Vit C (mg/dl)	3.6 ± 0.2	5.2 ± 0.2	5.3 ± 0.4	6.3 ± 0.2	p<0.05
Vit E (mg/L)	7.4 ± 0.6	8.2 ± 0.6	9.9 ± 0.7	11.2 ± 0.8	p<0.05
MDA (mmol/L)	2.8 ± 0.7	2.1 ± 0.5	1.3 ± 0.1	0.9 ± 0.1	p<0.05

Table 3. Comparison of oxidants and antioxidants profile in suspected cases of IEM.

Table 4. Correlations between oxidants and antioxidants in Suspected IEM cases.

Groups (MDA& TAC)	r-value	p value
High Risk	-0.015	p<0.05
Moderate Risk	-0.075	p<0.05
Low Risk	0.136	p<0.05
Controls	0.402	p<0.05

Note: *r*-value describes Negative and positive correlations between MDA and TAC activity.

groups, and positive correlations between low risk and control groups, suggested of oxidative stress. Our results clearly shows that severe oxidative stress in the IEM positive cases is compared with moderate and low risk category and increases significantly with the progression of disease, (i.e) the risk of IEM represents as high- risk positive patients for the specific diseases.

Discussion

The inborn errors of metabolism involve three inter-related and often overlapping factors as accumulation of toxic metabolite, deficiency of essential metabolites, and impaired functions of mitochondria. Oxidative stress often results has been proposed as a prime cause of the cellular injury that causes in the IEM. We quantified total antioxidant activity as early biochemical marker in various IEM to test their potential utility as clinically relevant marker. Likely SOD levels, Vit E, Vit C, MDA levels have been employed to measure oxidative and antioxidants stress in vivo using dried blood spots sampling as a secondary use of residual blood spots after it used for the neonatal screening.

In this study the biomarkers, taken as a whole, elevated markers of oxidative damage and depleted antioxidants defenses were seen in most of the IEM diagnostic group we assessed. In aminoaciduria's, the ketogenic acids have been showing to include oxidative models. Previous studies on Urea cycle disorders under treatment found elevated levels of oxidative damage were seen in subjects with endocrine metabolic disorders like CH and CAH. Decrease antioxidant activity was marked in organic acidurias[9-12]. The summary of the oxidative stress, where elevated levels in high risk patients compared with control groups. DBS-TAC can be used as a early biochemical marker for oxidative stress in high risk for IEM (Inborn errors of metabolism) which may result in reduced tissue damage by free radicals and help to monitor and optimize antioxidant therapy in such high risk neonates (Table 4).

The risk of IEM is inversely proportional to TAC in neonates. Under normal conditions, the ROS produced in the course of metabolism are contained by the natural antioxidant system which consists of a series of antioxidant enzymes such as superoxide dismutase (SOD), catalase, glutathione peroxidase as well as numerous endogenous and dietary antioxidant compounds that are capable of reacting with and inactivating ROS thereby protects the functional and structural molecules against ROS-mediated tissue damage. MDA is a three carbon, low molecular weight aldehyde that can be produced from free radical attack on polyunsaturated fatty acids of biological membranes [13,14].

The determination of MDA is used for monitoring lipid peroxidation in biological samples. TAC reflects the increased formation of ROS and lipid peroxidation with progressive symptoms to evalute. Vitamin E, a potent chain breaking lipid soluble antioxidant, reacts with lipid peroxyl radicals eventually terminating the peroxidation chain reaction and thereby reducing oxidative damage. Vitamin C represents the major water-soluble antioxidant in the human body. SOD is an endogenous antioxidant that catalyzes the disputation of the superoxide anion radical [15-19].There are conflicting reports in the values of antioxidant vitamins (vitamin E and vitamin C) and SOD enzyme activity among IEM high risk cases in various stages in the literature. In our study, vitamin E, vitamin C, SOD and TAC levels are decreased in high risk cases.

Conclusion

FRAP assay is presented as a novel method of assessing total antioxidant capacity and is considered as a useful indicator of the body's antioxidant status to counteract the oxidative damage due to ROS. The advantage of the FRAP assay is in being fast, easy to handle, with highly reproducible results. A significantly lower DBS TAC concentration in the IEM high risk patients compared to controls reflects a lower total antioxidant capacity. Thus, TAC may be useful as an early marker of oxidative stress to monitor and optimize antioxidant therapy as an adjunct in the management of IEM patients.

References

- 1. Seema K & Madhulika K. Newborn screening in india: Current perspectives. Indian Pediatr. 2010;47(3):219-4.
- 2. De Jesus VR, Zhang XK, Keutzer J, et al. Development and evaluation of quality control dried blood spot materials in

newborn screening for lysosomal storage disorders. Clin Chem. 2009;55(1):158-64.

- 3. Iris Benzie FF, Strain JJ. The ferric reducing ability of plasma (FRAP) as a measure of antioxidant power: the FRAP assay. Analytical Biochem. 1996;239(1):70-6.
- 4. Baker H, Frank O. Clinical vitaminology Academic Press. New York; 1969:169-73.
- Meltzer PS, Kallioniemi A, Trent JM. Chromosome alterations in human solid tumors. The genetic basis of human cancer. New York: McGraw-Hill. 2002:93-113.
- Omaye ST. Ascorbic acid analysis. II. Determination after derivatisation with 2.2. dinitrophenylhidrazine. Selected methods for determination of ascorbic acid in animal cells tissues and fluids. Methods in enzymology.1979;62:7-8.
- 7. Marklund SL. Spectrophotometric study of spontaneous disproportionation of superoxide anion radical and sensitive direct assay for superoxide dismutase. J Biol Chem.1976:251(23):7504-7.
- 8. Marklund SL. Direct assay with potassium superoxide. In Handbook of Methods for Oxygen Radical Research Edite
- 9. Dos Santos Mello M, Ribas GS, Wayhs CA, et al. Increased oxidative stress in patients with 3-hydroxy-3-methylglutaric aciduria. Mol Cell Biochem. 2015;402(1):149-55.
- 10. Guerreiro G, Mescka CP, Sitta A, et al. Urinary biomarkers of oxidative damage in Maple syrup urine disease: the L-carnitine role. Int J Dev Neurosci. 2015;42:10-4.
- 11. Mc Guire PJ, Parikh A, Diaz GA. Profiling of oxidative stress in patients with inborn errors of metabolism. Mol Genet Metab. 2009;98(1-2):173-80.

- 12. De Sousa MC, Vieira RB, Dos Santos DS, et al. Antioxidants and biomarkers of oxidative damage in the saliva of patients with Down's syndrome. Arch Oral Biol. 2015;60(4):600-5.
- Livnat G, Bentur L, Kuzmisnsky E, et al. Salivary profile and oxidative stress in children and adolescents with cystic fibrosis. J Oral Pathol Med. 2010;39(1):16-21
- Subramaniam P, Girish Babu K, Mohan Das L. Assessment of salivary total antioxidant levels and oral health status in children with Down syndrome. Spec Care Dentist. 2014;34(4):193-200.
- 15. Campos C, Guzman R, Lopez-Fernandez E, et al. Urinary uric acid and antioxidant capacity in children and adults with Down syndrome. Clin Biochem. 2010;43(3):228-33.
- 16. Ribas GS, Biancini GB, Mescka C, et al. Oxidative stress parameters in urine from patients with disorders of propionate metabolism: a beneficial effect of L:-carnitine supplementation. Cell Mol Neurobiol. 2012;32(1):77-82.
- 17. Cao G, Prior RL. Comparison of different analytical methods for assessing total antioxidant capacity of human serum. Clin Chem. 1998;44(6):1309-15.
- Bakhtiari A, Hajian-Tilaki K, Omidvar S, et al. Association of lipid peroxidation and antioxidant status with metabolic syndrome in Iranian healthy elderly women. Biomed Rep. 2017;7(4):331-6.
- 19. Ustundag Y, Huysal K, Kahvecioglu S, et al. Establishing reference values and evaluation of an in-house Ferric Reducing Antioxidant Power (FRAP) colorimetric assay in microplates. J Eur Res. 2016;2(2):126-31.